



CALAMBA WATER DISTRICT					
TITLE: OPERATION OF WISECLAVE® DIGITAL AUTOCLAVE					
DOC. NO. CWD-OPN-013		REVISION NO. 00		EFFECTIVE DATE: December 28, 2016	
Page 1 of 6					
		NAME		SIGNATURE	
AUTHOR		Ethel O. Paderes			
REVIEWED BY:		Engr. Joselito A. Gillera			
APPROVED BY:		Engr. Restituto B. Sumanga Sr.			
DOCUMENT HISTORY RECORD					
DCN	REV. NO.	DATE REVISED	AUTHOR	REASON FOR REVISION	
2016-12-033	00	N/A	Ethel O. Paderes	Initial Issue	

Important Note: This document and the information herein is the property of CALAMBA WATER DISTRICT, and issued in strict confidence. It shall not be reproduced, copied or given to a third party without express permission from CALAMBA WATER DISTRICT.														
Revision Locator						Master Copy Stamp		Copy Stamp		Copy Holder/ No. of Copies Issued				
Page #	1	2	3	4	5	6			1a			10		
Rev No.	00	00	00	00	00	00			1b			11		
Page #									2			12		
Rev No.									3			13		
Page #									4			14		
Rev No.									5			15		
Page #									6			16		
Rev No.									7			17		
Page #							8			18				
Rev No.							9			19				

CALAMBA WATER DISTRICT			
TITLE: OPERATION OF WISECLAVE® DIGITAL AUTOCLAVE			
DOC. NO. CWD-OPN-013	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 2 of 6

1.0 PURPOSE

- 1.1 The objective of this documented information is to provide standard instruction for the operation Wiseclave® Digital Autoclave
- 1.2 To determine and provide the resources needed to ensure valid and reliable results when monitoring or measuring is used to verify the conformity of products and services to requirements.
- 1.3 To ensure that the resources provided:
 - a) are suitable for the specific type of monitoring and measurement activities being undertaken;
 - b) are maintained to ensure continuing fitness for their purpose.
- 1.4 Retain appropriate documented information as evidence of fitness for purpose of the monitoring and measurement resources.

2.0 SCOPE


- 2.1 The scope applies to the safe operation of the autoclave by the authorized laboratory personnel.
- 2.2 Measurement traceability. When measurement traceability is a requirement, or is considered by CWD to be an essential part of providing confidence in the validity of measurement results, measuring equipment shall be:
 - a) calibrated or verified, or both, at specified intervals, or prior to use, against measurement standards traceable to international or national measurement standards; when no such standards exist, the basis used for calibration or verification shall be retained as documented information;
 - b) identified in order to determine their status;
 - c) safeguarded from adjustments, damage or deterioration that would invalidate the calibration status and subsequent measurement results.
- 2.3 To determine if the validity of previous measurement results has been adversely affected when measuring equipment is found to be unfit for its intended purpose, and shall take appropriate action as necessary.

3.0 RESPONSIBILITY

- 3.1 The Head of the Laboratory shall manage the administrative and technical operations of the laboratory.
- 3.2 The Laboratory Analyst must control and maintain the equipment.

4.0 DEFINITIONS


- 4.1 Sterilization - destruction or removal of all viable organisms from an object or from a particular environment
- 4.2 Autoclave - is a pressure cooker that sterilizes, or kills all microorganisms and their spores


Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder No. of Copies issued:	
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	1a	10
			1b	11
			2	12
			3	13
			4	14
			5	15
			6	16
			7	17
			8	18
			9	19

5.0 PROCESS FLOW STEPS

STEP	INPUT	PROCESS FLOW	OUTPUT	RESPONSIBLE
1. PREPARATION				
		START		
1	Media Preparation Logbook/ Equipment Monitoring Form	Fill the pressure vessel with distilled water up to the level of the heating coil		Laboratory Personnel
2		Ensure that the drain valve is closed prior to operation		Laboratory Personnel
3		Check the locking levers on the lid, it must be moved all the way into the locking ring		Laboratory Personnel
		END		
2. PROGRAM SETTING				
		START		
1		Press the MODE button		Laboratory Personnel
2		Change the position of the cursor by pressing the shift button		Laboratory Personnel
3		Change the minute setting using INC button		Laboratory Personnel
4		Press the MODE button to return to the main menu		Laboratory Personnel
5		Press Start / Stop button		Laboratory Personnel
		END		




<div>Master Copy Stamp:</div> <div><div><div>CALAMBA W.D.</div><div>MASTER COPY</div><div>12-28-16</div></div></div>	<div>Copy Stamp:</div>	<div>Important Note:</div> <div><div>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</div></div>	<div>Copy Holder/ No. of Copies issued:</div> <table><tr><td>1a</td><td>10</td></tr><tr><td>1b</td><td>11</td></tr><tr><td>2</td><td>12</td></tr><tr><td>3</td><td>13</td></tr><tr><td>4</td><td>14</td></tr><tr><td>5</td><td>15</td></tr><tr><td>6</td><td>16</td></tr><tr><td>7</td><td>17</td></tr><tr><td>8</td><td>18</td></tr><tr><td>9</td><td>19</td></tr></table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						

CALAMBA WATER DISTRICT			 <small>CALAMBA WATER DISTRICT</small>
TITLE: OPERATION OF WISECLAVE® DIGITAL AUTOCLAVE			
DOC. NO. CWD-OPN-013	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 4 of 6

STEP	INPUT	PROCESS FLOW	OUTPUT	RESPONSIBLE
1		<p>3. AFTER STERILIZATION</p> <p>4.</p> <pre> graph TD START([START]) --> OPEN[Open the autoclave after the pressure gauge has reach zero] OPEN --> END([END]) </pre>		Laboratory Personnel


6.0 PROCESS DETAILS:


6.1 Preparation





Fill the pressure vessel with distilled water up to the level of the heating coil	
Ensure that the drain valve is closed prior to operation	
Check the locking levers on the lid, it must be moved all the way into the locking ring	

6.2 Program Setting

Press the MODE button	
-----------------------	--

Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder No. of Copies issued:																				
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	<table><tr><td>1a</td><td>10</td></tr><tr><td>1b</td><td>11</td></tr><tr><td>2</td><td>12</td></tr><tr><td>3</td><td>13</td></tr><tr><td>4</td><td>14</td></tr><tr><td>5</td><td>15</td></tr><tr><td>6</td><td>16</td></tr><tr><td>7</td><td>17</td></tr><tr><td>8</td><td>18</td></tr><tr><td>9</td><td>19</td></tr></table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						

CALAMBA WATER DISTRICT			 CALAMBA WATER DISTRICT
TITLE: OPERATION OF WISECLAVE® DIGITAL AUTOCLAVE			
DOC. NO. CWD-OPN-013	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 5 of 6


Change the position of the cursor by pressing the shift button	
Change the minute setting using INC button	
Press the MODE button to return to the main menu	
Press Start / Stop button	


6.3 After sterilization

Open the autoclave after the pressure gauge has reach zero	
--	--

6.4 HEALTH AND SAFETY

- 6.4.1 Personal Protective Equipment (Apron, Goggles and Heat resistant gloves) must worn when removing the sterilized items
- 6.4.2 The handle should only be gripped at the experiment end and step back when opening the sterilizer
- 6.4.3 Autoclavable bags should be partially opened to prevent it from bursting.

Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder/ No. of Copies issued		
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	1a	10	
			1b	11	
			2	12	
			3	13	
			4	14	
			5	15	
			6	16	
			7	17	
			8	18	
			9	19	

CALAMBA WATER DISTRICT			 CALAMBA WATER DISTRICT
TITLE: OPERATION OF WISECLAVE® DIGITAL AUTOCLAVE			
DOC. NO. CWD-OPN-013	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 6 of 6

6.4.4 Do not sterilize pressure – sealed vessels

7.0 RECORDS RETENTION

7.1 Active Retention = indefinite retention period for current or active documents for both electronic and hardcopy/Controlled Copy.

7.2 Inactive/Archival Retention = shall be kept for active three (3) years or may request for an extension as deemed necessary (hardcopy); for electronic/soft file; it shall be kept in a separate folder named "Obsolete Copy".

8.0 REFERENCE DOCUMENTS

- 8.1 ISO 9001:2015 QMS Standard
- 8.2 QMS Manual
- 8.3 Statutory and Regulatory Requirements
- 8.4 WiseClave® Digital Autoclave Operation Manual
- 8.5 WHO Biosafety Manual 2004

9.0 ATTACHMENTS

- 9.1 Media Preparation Logbook
- 9.2 Equipment Monitoring Form


10.0 DISTRIBUTION LIST

Note 1: Select Relevant Recipient to Appear in below List.

COPY HOLDER NO.	DEPT/SEC./COPY HOLDER
1b	General Manager
8	Laboratory (Quality Control Division)

Note 2: Master Copy is in the custody of the Document Control Center.

- END -

Master Copy Stamp	Copy Stamp	Important Note	Copy Holder/ No. of Copies issued			
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	1a	10		
			1b	11		
			2	12		
			3	13		
			4	14		
			5	15		
			6	16		
			7	17		
			8	18		
			9	19		