





<b>CALAMBA WATER DISTRICT</b>					
<b>TITLE: Operation of Thermo Scientific Waterbath</b>					
DOCUMENT NO. CWD-OPN-001		REVISION NO. 00		EFFECTIVE DATE: December 28, 2016	
Page 1 of 6					
		NAME		SIGNATURE	
AUTHOR		Ethel O. Paderes			
REVIEWED BY:		Engr. Joselito A. Gillera			
APPROVED BY:		Engr. Restituto B. Sumanga Sr.			
<b>DOCUMENT HISTORY RECORD</b>					
DCN	REV. NO.	DATE REVISED	AUTHOR	REASON FOR REVISION	
2016-12-021	00	N/A	Ethel O. Paderes	Initial Issue	

**Important Note:**  
 This document and the information herein is the property of CALAMBA WATER DISTRICT and issued in strict confidence. It shall not be reproduced, copied or given to a third party without express permission from CALAMBA WATER DISTRICT.

Revision Locator							Master Copy Stamp:		Copy Stamp:		Copy Holder/ No. of Copies Issued			
Page #	1	2	3	4	5	6					1a	10		
Rev No.	00	00	00	00	00	00					1b	11		
Page #											2	12		
Rev No.											3	13		
Page #											4	14		
Rev No.											5	15		
Page #											6	16		
Rev No.											7	17		
Page #											8	18		
Rev No.											9	19		

<b>CALAMBA WATER DISTRICT</b>			
<b>TITLE: Operation of Thermo Scientific Waterbath</b>			
DOC. NO. CWD-OPN-001	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 2 of 6

## 1.0 PURPOSE

- 1.1 The objective of this procedure is to provide standard instruction for the operation Thermo Scientific Waterbath;
- 1.2 CWD shall implement production and service provision under controlled conditions. Controlled conditions shall include, as applicable:
- the availability of documented information that defines:
    - the characteristics of the products to be produced, the services to be provided, or the activities to be performed;
    - the results to be achieved;
    - the availability and use of suitable monitoring and measuring resources;
    - the implantation of monitoring and measurement activities at appropriate stages to verify that criteria for control of processes or outputs, and acceptance criteria for products and services, have been met;
    - the use of suitable infrastructure and environment for the operation of processes;
    - the appointment of competent persons, including any required qualification;
    - the validation, and periodic revalidation, of the ability to achieve planned results of the processes for production and service provision, where the resulting output cannot be verified by subsequent monitoring or measurement;
    - the implementation of actions to prevent human error;
    - the implementation of release, delivery and post-delivery activities.
  - To use suitable means to identify outputs when it is necessary to ensure the conformity of products and services.
  - To identify the status of outputs with respect to monitoring and measurement requirements throughout production and service provision.
  - To control the unique identification of the outputs when traceability is a requirement, and shall retain the documented information necessary to enable traceability.

## 2.0 SCOPE


- 2.1 The scope applies to the safe operation of the incubator by the authorized laboratory personnel.


## 3.0 RESPONSIBILITY

- 3.1 The Head of Laboratory shall manage the administrative and technical operations of the laboratory.
- 3.2 The Laboratory Analyst must control and maintain the equipment.

## 4.0 DEFINITION OF TERMS


- 4.1 Circulating Water Bath - water is thoroughly circulated throughout the bath resulting in a more uniform temperature.


Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder/ No. of Copies Issued:	
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	1a	10
			1b	11
			2	12
			3	13
			4	14
			5	15
			6	16
			7	17
			8	18
			9	19

<b>CALAMBA WATER DISTRICT</b>			
<b>TITLE: Operation of Thermo Scientific Waterbath</b>			
DOC. NO. CWD-OPN-001	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 3 of 6

#### 5.0 PROCESS FLOW STEPS




STEP	INPUT	PROCESS FLOW	OUTPUT	RESPONSIBLE
		START		
1	Operations Procedure	Add distilled water to bath		Laboratory Technician
2		Power up: Press the power switch located at the lower right hand corner of the bath		Laboratory Technician
3		Use the gable cover provided to conserve energy, reduce evaporation, and increase temperature control	Equipment Monitoring Form	Laboratory Technician
4		To set temperature control point		Laboratory Technician
5		Press SETPOINT. The display shows the current temperature set point		Laboratory Technician
6		Press UP and DOWN until the desired temperature set point displays		Laboratory Technician
7		Press ENTER to store the new value		Laboratory Technician
8		Display returns to actual temperature		Laboratory Technician
9		END		


Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder/ No. of Copies Issued:																				
<div><div>CALAMBA WD MASTER COPY</div><div>SIGNED: </div><div>DATE: 12-28-16</div></div>		<i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i>	<table><tr><td>1a</td><td>10</td></tr><tr><td>1b</td><td>11</td></tr><tr><td>2</td><td>12</td></tr><tr><td>3</td><td>13</td></tr><tr><td>4</td><td>14</td></tr><tr><td>5</td><td>15</td></tr><tr><td>6</td><td>16</td></tr><tr><td>7</td><td>17</td></tr><tr><td>8</td><td>18</td></tr><tr><td>9</td><td>19</td></tr></table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						

<b>CALAMBA WATER DISTRICT</b>			 ay Natinong
<b>TITLE: Operation of Thermo Scientific Waterbath</b>			
DOC. NO. CWD-OPN-001	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 4 of 6

## 6.0 PROCESS DETAILS

### 6.1 Operation

<p>Add distilled water to bath</p> <ul style="list-style-type: none"> <li>Water Level Requirements</li> </ul> <p>Minimum level – 1 ½ inches above the pump outlet</p>	
<p>Power up: Press the power switch located at the lower right hand corner of the bath. The unit is fully operational upon application of power and the display will read actual bath temperature</p>	
<p>Use the gable cover provided to conserve energy, reduce evaporation, and increase temperature control accuracy.</p>	


Master Copy Stamp:	Copy Stamp:	Important Note	Copy Holder No. of Copies Issued	
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	1a	10
			1b	11
			2	12
			3	13
			4	14
			5	15
			6	16
			7	17
			8	18
			9	19

<b>CALAMBA WATER DISTRICT</b>			
<b>TITLE: Operation of Thermo Scientific Waterbath</b>			
DOC. NO. CWD-OPN-001	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 5 of 6

## 6.2 To set temperature control point

- a. Press SETPOINT. The display shows the current temperature set point.
- b. Press UP and DOWN until the desired temperature set point displays.
- c. Press ENTER to store the new value.
- d. Display returns to actual temperature

\*If enter is not pressed, new value will not be set and display will return to actual temperature



## 6.3 HEALTH AND SAFETY

- 6.3.1 Personal Protective Equipment (Apron, Goggles and Heat resistant gloves) must be worn when placing and removing samples for incubation.
- 6.3.2 Care should be exercised when removing the back cover. The cover should be placed upside down when removed to maintain a dry work area.
- 6.3.3 Do not touch the metal surface of the gable cover when operating at temperatures greater than 60°C.
- 6.3.4 Avoid spilling harsh chemicals onto the bath, as corrosion of the stainless steel may result.

## 7.0 RECORDS RETENTION


- 7.1 Active Retention – indefinite retention period for current or active documents for both electronic and hardcopy Master Copy.
- 7.2 Inactive/Archival Retention – shall be kept for active three (3) years or may request for an extension as deemed necessary (hardcopy); for electronic/soft file; it shall be kept in a separate folder named "Obsolete Master Copy/Original".

## 8.0 REFERENCE


- 8.1 ISO 9001:2015 QMS Standard
- 8.2 Quality Manual
- 8.3 Thermo Scientific Waterbath Manual

## 9.0 ATTACHMENTS

- 9.1 Equipment Monitoring Form

Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder No. of Copies Issued	
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	1a	10
			1b	11
			2	12
			3	13
			4	14
			5	15
			6	16
			7	17
			8	18
			9	19



<b>CALAMBA WATER DISTRICT</b>			
<b>TITLE: Operation of Thermo Scientific Waterbath</b>			
DOC. NO. CWD-OPN-001	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 6 of 6


#### 10.0 DISTRIBUTION LIST

Note 1: Select Relevant Recipient to Appear in below List.

COPY HOLDER NO.	DEPT/SEC./COPY HOLDER
1b	General Manager
8	Laboratory (Quality Control Division)

Note 2: Master Copy is in the custody of the Document Control Center.

- END

Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder/ No. of Copies Issued:			
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	1a		10	
			1b		11	
			2		12	
			3		13	
			4		14	
			5		15	
			6		16	
			7		17	
			8		18	
			9		19	