


<b>CALAMBA WATER DISTRICT</b>					
<b>TITLE: OPERATION OF REXMED AUTOMATIC AUTOCLAVE</b>					
DOCUMENT NO. CWD-OPN-003		REVISION NO. 00		EFFECTIVE DATE: December 28, 2016	
Page 1 of 6					
		NAME		SIGNATURE	
AUTHOR		Ethel O. Paderes			
REVIEWED BY:		Engr. Joselito A. Gillera			
APPROVED BY:		Engr. Restituto B. Sumanga SR.			
<b>DOCUMENT HISTORY RECORD</b> <span style="float: right;">Form No.</span>					
DCN	REV. NO.	DATE REVISED	AUTHOR	REASON FOR REVISION	
2016-12-023	00	N/A	Ethel O. Paderes	Initial Issue	

**Important Note:**  
 This document and the information herein is the property of CALAMBA WATER DISTRICT and issued in strict confidence. It shall not be reproduced, copied or given to a third party without express permission from CALAMBA WATER DISTRICT.

Revision Locator							Master Copy Stamp:		Copy Stamp:		Copy Holder/ No. of Copies Issued				
Page #	1	2	3	4	5	6				1a			10		
Rev No.	00	00	00	00	00	00				1b			11		
Page #										2			12		
Rev No.										3			13		
Page #										4			14		
Rev No.										5			15		
Page #										6			16		
Rev No.										7			17		
Page #										8			18		
Rev No.										9			19		

<b>CALAMBA WATER DISTRICT</b>			
<b>TITLE: OPERATION OF REXMED AUTOMATIC AUTOCLAVE</b>			
DOC. NO. CWD-OPN-003	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 2 of 6

## 1.0 PURPOSE

- 1.1 The objective of this documented information is to provide standard instruction for the operation of Rexmed Automatic Autoclave LS – 2D.
- 1.2 CWD shall implement production and service provision under controlled conditions. Controlled conditions shall include, as applicable:
- the availability of documented information that defines:
    - the characteristics of the products to be produced, the services to be provided, or the activities to be performed;
    - the results to be achieved;
  - the availability and use of suitable monitoring and measuring resources;
  - the implantation of monitoring and measurement activities at appropriate stages to verify that criteria for control of processes or outputs, and acceptance criteria for products and services, have been met;
  - the use of suitable infrastructure and environment for the operation of processes;
  - the appointment of competent persons, including any required qualification;
  - the validation, and periodic revalidation, of the ability to achieve planned results of the processes for production and service provision, where the resulting output cannot be verified by subsequent monitoring or measurement;
  - the implementation of actions to prevent human error;
  - the implementation of release, delivery and post-delivery activities.
- 1.3 To use suitable means to identify outputs when it is necessary to ensure the conformity of products and services.
- 1.4 To identify the status of outputs with respect to monitoring and measurement requirements throughout production and service provision.
- 1.5 To control the unique identification of the outputs when traceability is a requirement, and shall retain the documented information necessary to enable traceability.

## 2.0 SCOPE


- 2.1 The scope applies to the safe operation of the autoclave by the authorized laboratory personnel.

## 3.0 RESPONSIBILITY


- 3.1 The Head of Laboratory shall manage the administrative and technical operations of the laboratory.
- 3.2 The Laboratory Analyst must control and maintain the equipment.

## 4.0 DEFINITION OF TERMS

- 4.1 Sterilization – the removal of all microorganism and pathogens from an object or surface by treating it with chemicals by subjecting it to high heat or radiation.

Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder No. of Copies Issued		
		<i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i>	1a	10	
			1b	11	
			2	12	
			3	13	
			4	14	
			5	15	
			6	16	
			7	17	
			8	18	
			9	19	

**5.0 PROCESS FLOW STEPS**



Master Copy Stamp: <div style="border: 2px solid green; padding: 10px; text-align: center; margin-top: 10px;"> <b>CALAMBAWD</b>  <b>MASTER</b>  <b>COPY</b>                      SIGNED:                       DATE: 12-28-16                 </div>	Copy Stamp:	Important Note: <p style="margin-top: 20px;"><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	Copy Holder/ No. of Copies Issued: <table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <tr> <td>1a</td> <td>10</td> </tr> <tr> <td>1b</td> <td>11</td> </tr> <tr> <td>2</td> <td>12</td> </tr> <tr> <td>3</td> <td>13</td> </tr> <tr> <td>4</td> <td>14</td> </tr> <tr> <td>5</td> <td>15</td> </tr> <tr> <td>6</td> <td>16</td> </tr> <tr> <td>7</td> <td>17</td> </tr> <tr> <td>8</td> <td>18</td> </tr> <tr> <td>9</td> <td>19</td> </tr> </table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						

<b>CALAMBA WATER DISTRICT</b> <b>TITLE: OPERATION OF REXMED AUTOMATIC AUTOCLAVE</b>			
DOC. NO. CWD-OPN-003	REVISION NO. 00	EFFECTIVE DATE. December 28, 2016	

Page 4 of 6

## 6.0 PROCESS DETAILS

### 6.1 Preparation

Open the lid and remove the bottom plate. Fill the pressure vessel with distilled water about 1.5 liters, the level must cover the top of heating element.



Return the bottom plate.




Place the item to be sterilized inside the chamber. Close the lid securely and tightly.



### 6.2 Program Setting

Put the selector switch to "Ster" position





Master Copy Stamp	Copy Stamp	Important Note	Copy Holder/ No. of Copies Issued																				
<div><div>CALAMBAWD MASTER COPY</div><div>SIGNED: </div><div>DATE: 12-28-16</div></div>		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	<table><tr><td>1a</td><td>10</td></tr><tr><td>1b</td><td>11</td></tr><tr><td>2</td><td>12</td></tr><tr><td>3</td><td>13</td></tr><tr><td>4</td><td>14</td></tr><tr><td>5</td><td>15</td></tr><tr><td>6</td><td>16</td></tr><tr><td>7</td><td>17</td></tr><tr><td>8</td><td>18</td></tr><tr><td>9</td><td>19</td></tr></table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						


<b>CALAMBA WATER DISTRICT</b> <b>TITLE: OPERATION OF REXMED AUTOMATIC AUTOCLAVE</b>			
DOC. NO. CWD-OPN-003	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	

Page 5 of 6

Set timer for sterilization 15 mins – media 30 mins - glassware		
Put master switch "On", green pilot lamp glows		
The red pilot indicates timer countdown		

### 6.3 After sterilization

The buzzer will sound after the sterilization cycle		
Turn Master Switch to "off"		

Master Copy Stamp:	Copy Stamp	Important Note:	Copy Holder's No. of Copies Issued
		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	1a
			1b
			2
			3
			4
			5
			6
			7
			8
			9
			10
			11
			12
			13
			14
			15
			16
			17
			18
			19

<b>CALAMBA WATER DISTRICT</b>			 CAL Water District
<b>TITLE: OPERATION OF REXMED AUTOMATIC AUTOCLAVE</b>			
DOC. NO. CWD-OPN-003	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 6 of 6

Open the exhaust knob valve for steam exhaustion until 0 pressure



#### 6.4 HEALTH AND SAFETY

- 6.4.1 Personal Protective Equipment (Apron, Goggles and Heat resistant gloves) must worn when removing the sterilized items
- 6.4.2 The handle should only be gripped at the experiment end and step back when opening the sterilizer
- 6.4.3 Autoclavable bags should be partially opened to prevent it from bursting.
- 6.4.4 Do not sterilize pressure – sealed vessels the bath, as corrosion of the stainless steel may result.

#### 7.0 RECORDS RETENTION

- 7.1 Active Retention – indefinite retention period for current or active documents for both electronic and hardcopy Master Copy.
- 7.2 Inactive/Archival Retention – shall be kept for active three (3) years or may request for an extension as deemed necessary (hardcopy); for electronic/soft file; it shall be kept in a separate folder named "Obsolete Master Copy/Original".

#### 8.0 REFERENCE

- 8.1 ISO 9001:2015 QMS Standard
- 8.2 Rexmed Automatic Autoclave Operation Manual
- 8.3 WHO Biosafety Manual 2004

#### 9.0 ATTACHMENTS

- 9.1 Equipment Monitoring Form


#### 10.0 DISTRIBUTION LIST

**Note 1: Select Relevant Recipient to Appear in below List.**

COPY HOLDER NO.	DEPT/SEC./COPY HOLDER
1b	General Manager
8	Laboratory (Quality Control Division)

**Note 2: Master Copy is in the custody of the Document Control Center.**

- END

Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder/ No. of Copies Issued:																				
<div><div>CALAMBAWD MASTER COPY</div><div>SIGNED: </div><div>DATE: 12-28-16</div></div>		<p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	<table><tr><td>1a</td><td>10</td></tr><tr><td>1b</td><td>11</td></tr><tr><td>2</td><td>12</td></tr><tr><td>3</td><td>13</td></tr><tr><td>4</td><td>14</td></tr><tr><td>5</td><td>15</td></tr><tr><td>6</td><td>16</td></tr><tr><td>7</td><td>17</td></tr><tr><td>8</td><td>18</td></tr><tr><td>9</td><td>19</td></tr></table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						