





CALAMBA WATER DISTRICT					
TITLE: OPERATION AND CALIBRATION OF INOLAB PH LEVEL 1					
DOCUMENT NO. CWD-OPN-006		REVISION NO. 00		EFFECTIVE DATE: December 28, 2016	
Page 1 of 6					
		NAME		SIGNATURE	
AUTHOR		Ethel O. Paderes			
REVIEWED BY:		Engr. Joselito A. Gillera			
APPROVED BY:		Engr. Restituto B. Sumanga Sr.			
DOCUMENT HISTORY RECORD Form No.					
DCN	REV. NO.	DATE REVISED	AUTHOR	REASON FOR REVISION	
2016-12-026	00	N/A	Ethel O. Paderes	Initial Issue	

Important Note:
This document and the information herein is the property of CALAMBA WATER DISTRICT and issued in strict confidence. It shall not be reproduced, copied or given to a third party without express permission from CALAMBA WATER DISTRICT.

Revision Locator								Master Copy Stamp:		Copy Stamp:		Copy Holder/ No. of Copies Issued			
Page #	1	2	3	4	5	6						1a		10	
Rev No.	00	00	00	00	00	00						1b		11	
Page #												2		12	
Rev No.												3		13	
Page #												4		14	
Rev No.												5		15	
Page #												6		16	
Rev No.												7		17	
Page #								8		18					
Rev No.								9		19					

CALAMBA WATER DISTRICT TITLE: OPERATION AND CALIBRATION OF INOLAB pH LEVEL 1			
DOC. NO. CWD-OPN-006	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	

Page 2 of 6

1.0 PURPOSE

- 1.1 The objective of this documented information is to provide standard instruction for the operation and calibration of inoLab Level 1.
- 1.2 To determine and provide the resources needed to ensure valid and reliable results when monitoring or measuring is used to verify the conformity of products and services to requirements.
- 1.3 To ensure that the resources provided:
 - a) are suitable for the specific type of monitoring and measurement activities being undertaken;
 - b) are maintained to ensure continuing fitness for their purpose.
- 1.4 retain appropriate documented information as evidence of fitness for purpose of the monitoring and measurement resources.

2.0 SCOPE


- 2.1 The scope applies to the safe operation of the pH meter by the authorized laboratory personnel.
- 2.2 Measurement traceability. When measurement traceability is a requirement, or is considered by CWD to be an essential part of providing confidence in the validity of measurement results, measuring equipment shall be:
 - a) calibrated or verified, or both, at specified intervals, or prior to use, against measurement standards traceable to international or national measurement standards; when no such standards exist, the basis used for calibration or verification shall be retained as documented information;
 - b) identified in order to determine their status;
 - c) safeguarded from adjustments, damage or deterioration that would invalidate the calibration status and subsequent measurement results.
- 2.3 To determine if the validity of previous measurement results has been adversely affected when measuring equipment is found to be unfit for its intended purpose, and shall take appropriate action as necessary.


3.0 RESPONSIBILITY

- 3.1 The Head of Laboratory shall manage the administrative and technical operations of the laboratory.
- 3.2 The Laboratory Analyst must control and maintain the equipment.

4.0 DEFINITION OF TERMS


- 4.1 Calibration – is the comparison of measurement values delivered by a device under test with those of a calibration standard of known accuracy.
- 4.2 pH - a numeric scale used to specify the acidity or basicity of an aqueous solution.


Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder/ No. of Copies Issued:																				
		<i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i>	<table><tr><td>1a</td><td>10</td></tr><tr><td>1b</td><td>11</td></tr><tr><td>2</td><td>12</td></tr><tr><td>3</td><td>13</td></tr><tr><td>4</td><td>14</td></tr><tr><td>5</td><td>15</td></tr><tr><td>6</td><td>16</td></tr><tr><td>7</td><td>17</td></tr><tr><td>8</td><td>18</td></tr><tr><td>9</td><td>19</td></tr></table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						




CALAMBA WATER DISTRICT TITLE: OPERATION AND CALIBRATION OF INOLAB pH LEVEL 1			
DOC. NO. CWD-OPN-006	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	


Page 3 of 6


5.0 PROCESS FLOW STEPS

STEP	INPUT	PROCESS FLOW	OUTPUT	RESPONSIBLE
				
1	Inolab Operation Procedure	Place the instrument on a flat surface and protect it against intense light and heat		Laboratory Analyst
2		Press the key		Laboratory Analyst
3		Connect the pH electrode to the instrument. The measuring instrument is ready for operation		Laboratory Analyst
4		Calibrating (pH) ConCal		Laboratory Analyst
5		With select the pH mode and with select the ConCal procedure.		Laboratory Analyst
6		Submerge the pH electrode in the pH 7.0 ± 0.5 buffer solution.		Laboratory Analyst
7		Press the key. The measured pH value appears on the display.		Laboratory Analyst
8		Press the key. The value of the asymmetry (mV) and the sensor symbol appear on the display.		Laboratory Analyst
9		Press the key. SLO(pe) appears on the display.		Laboratory Analyst
		Thoroughly rinse the electrode with distilled water.		Laboratory Analyst
		Submerge the pH electrode in the second buffer solution.		Laboratory Analyst
		A		

Master Copy Stamp:	Copy Stamp:	Important Note:	Copy Holder/ No. of Copies Issued:																				
		<i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i>	<table><tr><td>1a</td><td>10</td></tr><tr><td>1b</td><td>11</td></tr><tr><td>2</td><td>12</td></tr><tr><td>3</td><td>13</td></tr><tr><td>4</td><td>14</td></tr><tr><td>5</td><td>15</td></tr><tr><td>6</td><td>16</td></tr><tr><td>7</td><td>17</td></tr><tr><td>8</td><td>18</td></tr><tr><td>9</td><td>19</td></tr></table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						

CALAMBA WATER DISTRICT				
TITLE: OPERATION AND CALIBRATION OF INOLAB pH LEVEL 1				
DOC. NO. CWD-OPN-006		REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	
		Page 4 of 6		
STEP	INPUT	PROCESS FLOW	OUTPUT	RESPONSIBLE
		 <pre> graph TD A((A)) --> S1[Press the key. The second measured pH value appears on the display.] S1 --> S2[Press the key. The value of the slope (mV/pH) appears on the display.] S2 --> S3[Press the key. The value of the asymmetry (mV) appears on the display again.] S3 --> S4[To return to the measuring mode: Press the key.] S4 --> S5[Connect the electrode to the pH meter] S5 --> S6[Switch on the meter by pressing the key] S6 --> S7[Calibrate the measuring system.] S7 -.-> E1[Equipment Monitoring Form] S7 --> S8[Submerge the electrode to the test sample.] S8 --> S9[Select the measuring mode by pressing the MODE key] S9 --> S10[Wipe the outside of the measuring instrument with a damp, lint free cloth. Disinfect the housing with 70%] S10 --> END((END)) </pre>		
1		Press the key. The second measured pH value appears on the display.		Laboratory Analyst
2		Press the key. The value of the slope (mV/pH) appears on the display.		Laboratory Analyst
3		Press the key. The value of the asymmetry (mV) appears on the display again.		Laboratory Analyst
4		To return to the measuring mode: Press the key.		Laboratory Analyst
5		Connect the electrode to the pH meter		Laboratory Analyst
6		Switch on the meter by pressing the key		Laboratory Analyst
7		Calibrate the measuring system.		Laboratory Analyst
8		Submerge the electrode to the test sample.		Laboratory Analyst
9		Select the measuring mode by pressing the MODE key		Laboratory Analyst
		Wipe the outside of the measuring instrument with a damp, lint free cloth. Disinfect the housing with 70%		Laboratory Analyst
		END		


Master Copy Stamp: 	Copy Stamp: 	Important Note: <p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	Copy Holder/ No. of Copies Issued: <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>1a</td><td>10</td></tr> <tr><td>1b</td><td>11</td></tr> <tr><td>2</td><td>12</td></tr> <tr><td>3</td><td>13</td></tr> <tr><td>4</td><td>14</td></tr> <tr><td>5</td><td>15</td></tr> <tr><td>6</td><td>16</td></tr> <tr><td>7</td><td>17</td></tr> <tr><td>8</td><td>18</td></tr> <tr><td>9</td><td>19</td></tr> </table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						

CALAMBA WATER DISTRICT TITLE: OPERATION AND CALIBRATION OF INOLAB pH LEVEL 1			
DOC. NO. CWD-OPN-006	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	

Page 5 of 6

6.0 PROCESS DETAILS

6.1 Installation

- 6.1.1 Place the instrument on a flat surface and protect it against intense light and heat
- 6.1.2 Press  the key.

The Display test appears briefly on the display. The instrument then switches automatically to the previously selected measuring mode.

- 6.1.3 Connect the pH electrode to the instrument. The measuring instrument is ready for operation.

6.2. Calibrating (pH) ConCal


Use two buffer solutions for this procedure:

- pH 7.0 ± 0.5
- any other buffer solution

- 6.2.1 With  select the pH mode and  with select the ConCal procedure.

- 6.2.2 Submerge the pH electrode in the pH 7.0 ± 0.5 buffer solution.

- 6.2.3 Press the  key. The measured pH value appears on the display.


- 6.2.4 Press the  key. The value of the asymmetry (mV) and the sensor symbol appear on the display.

- 6.2.5 Press the  key. SLO(pe) appears on the display.

- 6.2.6 Thoroughly rinse the electrode with distilled water.

- 6.2.7 Submerge the pH electrode in the second buffer solution.

- 6.2.8 Press the  key. The second measured pH value appears on the display.

- 6.2.9 Press the  key. The value of the slope (mV/pH) appears on the display.

The sensor symbol shows the evaluation of the electrode after the two-point calibration.

- 7.0 Press the  key. The value of the asymmetry (mV) appears on the display again.

- 7.1 To return to the measuring mode: Press the  key.


6.3 Measuring pH


- 6.3.1 Connect the electrode to the pH meter.

- 6.3.2 Switch on the meter by pressing .

- 6.3.3 Calibrate the measuring system.

- 6.3.4 Submerge the electrode to the test sample.

- 6.3.5 Select the measuring mode by pressing . The display shows the pH value or the Redox voltage (mV) of the sample.

Master Copy Stamp: 	Copy Stamp:	Important Note: <i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i>	Copy Holder/ No. of Copies Issued: <table><tr><td>1a</td><td>10</td></tr><tr><td>1b</td><td>11</td></tr><tr><td>2</td><td>12</td></tr><tr><td>3</td><td>13</td></tr><tr><td>4</td><td>14</td></tr><tr><td>5</td><td>15</td></tr><tr><td>6</td><td>16</td></tr><tr><td>7</td><td>17</td></tr><tr><td>8</td><td>18</td></tr><tr><td>9</td><td>19</td></tr></table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						

CALAMBA WATER DISTRICT			 CWD Logo Banner
TITLE: OPERATION AND CALIBRATION OF INOLAB pH LEVEL 1			
DOC. NO. CWD-OPN-006	REVISION NO. 00	EFFECTIVE DATE: December 28, 2016	Page 6 of 6

6.4 Cleaning and Disinfection

6.4.1 Wipe the outside of the measuring instrument with a damp, lint free cloth. Disinfect the housing with 70% alcohol.

7.0 RECORDS RETENTION

7.1 Active Retention – indefinite retention period for current or active documents for both electronic and hardcopy Master Copy.

7.2 Inactive/Archival Retention – shall be kept for active three (3) years or may request for an extension as deemed necessary (hardcopy); for electronic/soft file; it shall be kept in a separate folder named "Obsolete Master Copy/Original".

8.0 REFERENCE

8.1 ISO 9001:2015 QMS Standard

8.2 Quality Manual (Optional)

8.3 Thermo Scientific Waterbath Manual

9.0 ATTACHMENTS

9.1 Equipment Monitoring Form


10.0 DISTRIBUTION LIST

Note 1: Select Relevant Recipient to Appear in below List.

COPY HOLDER NO.	DEPT/SEC./COPY HOLDER
1b	General Manager
8	Laboratory (Quality Control Division)

Note 2: Master Copy is in the custody of the Document Control Center.

- END

Master Copy Stamp: 	Copy Stamp: 	Important Note: <p><i>This documented information is not to be reproduced in any form without permission; and shall not be discarded unless superseded by a revised issue.</i></p>	Copy Holder/ No. of Copies Issued: <table border="1"> <tr><td>1a</td><td>10</td></tr> <tr><td>1b</td><td>11</td></tr> <tr><td>2</td><td>12</td></tr> <tr><td>3</td><td>13</td></tr> <tr><td>4</td><td>14</td></tr> <tr><td>5</td><td>15</td></tr> <tr><td>6</td><td>16</td></tr> <tr><td>7</td><td>17</td></tr> <tr><td>8</td><td>18</td></tr> <tr><td>9</td><td>19</td></tr> </table>	1a	10	1b	11	2	12	3	13	4	14	5	15	6	16	7	17	8	18	9	19
1a	10																						
1b	11																						
2	12																						
3	13																						
4	14																						
5	15																						
6	16																						
7	17																						
8	18																						
9	19																						